

# OAPL ANKLE WALKER

A Guide for Patients

 oapl.™



# INTRODUCTION

You have been fitted with an OAPL Ankle Walker to control the alignment and reduce the movement of the bones and joints in your lower leg and foot.

The OAPL Ankle Walker is available in 2 styles:

1. Standard Ankle Walker with a Foam Liner
2. Air Ankle Walker with an Inflationable Liner

Patient's are usually fitted with an ankle walker due to one of the following:

- Post operatively due to foot or ankle surgery.
- To restrict movement following soft tissue injuries of the ankle.
- Fractures and/or Achilles tendon repair or injury.

**Consulting Orthotist**

---

**Clinic Attended**

---

**Date Attended**

---



# APPLICATION

Your Orthotist will fit the walker to your injured limb. It will be set up in the ideal position for your injury.

## Step 1

Leave the liner attached to the plastic uprights. Put your leg into the walker ensuring that your heel is right to the back and you can feel it contacting with the bottom.



## Step 2

Close the liner firmly around your calf and across the top of your foot.



## Step 3

Fasten the straps.

Begin with the strap over the middle of your foot closest to the ankle.

Feed the strap through the "D" ring and pull it firmly. Use the pads provided under this to prevent rubbing.



## Step 4

Next, fasten the strap closest to the ankle.

Followed by the two remaining straps on the leg.

Ensure that all straps are secured firmly. They should not cause any discomfort but need to be firm in order to immobilise your foot and ankle.

If you are unsure about how firm they need to be, get your Orthotist to mark the straps so that there is no doubt.



# APPLICATION

## Step 5

Finally, fasten the last strap across the foot.

If you have an inflatable liner model turn the valve in a clockwise direction and use the blue button to inflate the liner until a snug fit is achieved.

## Step 6

To remove the walker:

1. Deflate liner (turn valve anti-clockwise direction).
2. Undo all the straps in any order making sure that you stick the velcro back on itself. This prevents it from getting tangled and sticking to everything else (i.e. clothes).
3. Open up the liner down the front of your leg.
4. Open up the liner over the foot.
5. Lift your leg out leaving liner attached to the plastic base and uprights.



**Check with your specialist or orthotist as to how much weight you can take on the affected side.**

# WEIGHT BEARING INSTRUCTIONS

You may find your walker a little bit difficult to walk in at first until you get used to it. If you have been using crutches or a stick you may want to continue to use these for a few days until you have adjusted.

The heel height of your walker is about the same as a pair of running shoes. It is therefore recommended that you wear a runner or equivalent on the other foot to try and keep you as even and balanced as possible.

It is suggested that you wear the walker with a sock so that the sock rather than the liner of the walker absorbs any perspiration. This will also keep your toes warm.

Your specialist/doctor will advise as to the length of time that you will be required to wear the ankle walker.

You should have a review appointment with them in several weeks time to assess your progress.

# CONTACT DETAILS

Email: [info@oapl.com.au](mailto:info@oapl.com.au)

## **Fitzroy Clinic**

93-97 Webb Street  
Fitzroy VIC 3065  
T: 9419 2499  
F: 9416 3543

## **Footscray Clinic**

55 Pickett Street  
Footscray VIC 3011  
T: 9687 4589  
F: 9687 3961

## **Clayton Clinic**

281 Clayton Road  
Clayton VIC 3168  
T: 9562 9422  
F: 9562 9605

## **Box Hill Clinic**

Epworth Eastern  
Suite 5A, Level 2  
1 Arnold Street  
Box Hill VIC 3128  
T: 9897 1374  
F: 9783 6944

## **Ringwood Clinic**

86 Mt Dandenong Road  
Ringwood East VIC 3135  
T: 9879 4299  
F: 9870 9411

## **Brighton Clinic**

Suite 4, 214 Bay Street  
Brighton VIC 3186  
T: 9596 6895  
F: 9824 8205

## **Frankston Clinic**

346 Nepean Highway  
Frankston VIC 3199  
T: 9783 3866  
F: 9783 6944

## **Richmond Clinic**

Epworth Centre  
Suite 5, Level 7  
32 Erin Street  
Richmond VIC 3121  
T: 9421 6226  
F: 9426 4321

## **Bendigo Clinic**

401 - 405 High Street  
Golden Square VIC 3555  
T: 5441 4333  
F: 5441 6555

## **Sydney Clinic**

Macquarie University Clinic 307  
Level 3, 2 Technology Place  
Macquarie University  
NSW 2109  
T: 02 9812 3854

29 South Corporate Ave  
Rowville Victoria 3178

T: 1300 TO OAPL  
or 1300 86 6275

F: 1300 OAPL FX  
or 1300 6275 39

[info@oapl.com.au](mailto:info@oapl.com.au)

ACN 006 520 314  
ABN 30 006 520 314

[www.oapl.com.au](http://www.oapl.com.au)

RPT0018  
Version 6.13

